

## Volumetric capnography option

(159136 CO<sub>2</sub> preparation kit, 281718 HAMILTON Capnostat 5 CO<sub>2</sub> sensor)

Capnography has become an important tool to evaluate the adequacy of ventilation, as the obvious goal of ventilation is to remove the volume of CO<sub>2</sub> produced by the body's metabolic processes. Conventionally, minute ventilation (MinVol) is used as a measure of overall ventilation. However, minute ventilation does not tell how much CO<sub>2</sub>-removing volume reaches the alveoli, because it includes not only ventilation of the lungs but also ventilation that is wasted in the airways. Thus, minute ventilation does not conclusively indicate the actual alveolar reach. Instead, alveolar ventilation (V'alv) should be measured\*.

The HAMILTON-G5 is the only ventilator that provides the option to measure alveolar ventilation plus other measures of volumetric capnography. With the unique combination of fast capnography and accurate spirometry, the HAMILTON-G5 accurately determines a number of other parameters:


- The CO<sub>2</sub> elimination (V'CO<sub>2</sub>) measurement permits assessment of metabolic rate (e.g., it is high with sepsis, tetanus, etc.) and treatment progress
- The end-tidal CO<sub>2</sub> (PetCO<sub>2</sub>) measurement permits assessment of arterial CO<sub>2</sub> partial pressure. (Note that it is inaccurate in pulmonary embolism)
- The airway dead space (VDaw) measurement permits assessment of actual alveolar ventilation (as opposed to minute ventilation)
- The capnogram shape (slopeCO<sub>2</sub>) permits assessment of COPD, asthma, inefficient ventilation, gastric intubation, and pulmonary embolism
- The physiological dead space fraction (VD/VT) permits assessment of risk [Nuckton 2002]

### Method

Alveolar ventilation is defined as minute ventilation minus dead space ventilation. Dead space is measured using the volumetric capnogram, which is divided into phase I (no CO<sub>2</sub> present), phase II (rapid rise in CO<sub>2</sub>), and phase III (alveolar plateau). Phases I and II represent airway dead space.

To obtain the volume of CO<sub>2</sub> exhaled in a single breath, the conventional CO<sub>2</sub> versus time curve and the spirogram are combined into the CO<sub>2</sub> versus volume diagram shown in the figure below. This requires the simultaneous measurement of flow and CO<sub>2</sub> at the airway opening.

The slope of phase III (slopeCO<sub>2</sub>) is a very sensitive lung function indicator. Like other parameters, it is calculated breath-by-breath on the HAMILTON-G5. A steep slope is seen in COPD patients, while a flat plateau is seen in post-operative patients. The parameter slopeCO<sub>2</sub> can be trended and effects of treatment thus be easily followed.



### Description of volumetric capnogram:

- Phase I: Pure airway dead space, from point of measurement of CO<sub>2</sub> toward the lungs
- Phase II: Weighted average of alveolar gas from different lung spaces, at the sensor location; measurement is VDaw
- Phase III: Alveolar plateau; measurement is slopeCO<sub>2</sub> together with end-tidal CO<sub>2</sub>.

\* For example, in a patient with airway dead space of 150 ml (68 kg IBW male), a tidal volume of 100 ml at a rate of 60 b/min yields the same minute ventilation as a tidal volume of 500 ml at a rate of 12 b/min, yet it has no real benefit to the patient since only dead space ventilation occurs.

## Ordering information

- 159136 CO<sub>2</sub> preparation kit (without sensor)  
 281718 HAMILTON CAPNOSTAT 5 CO<sub>2</sub> sensor  
 281719 CO<sub>2</sub> sensor airway adapter, adult/pediatric, single-use, box of 10  
 281720 CO<sub>2</sub> sensor airway adapter, infant/pediatric, single-use, box of 10  
 281721 CO<sub>2</sub> sensor airway adapter, adult/pediatric, reusable, box of 1  
 281722 CO<sub>2</sub> sensor airway adapter, infant/pediatric, reusable, box of 1  
 281803 15mm Male/Female adapter for infant flow sensor, package of 25

## Specifications for the CO<sub>2</sub> option

Parameter	Range	Resolution
FetCO <sub>2</sub>	0 to 19.7%	0.1%
PetCO <sub>2</sub>	0 to 150 mmHg	1 mmHg
slopeCO <sub>2</sub>	0 to 9.99 %CO <sub>2</sub> /l	0.01 %CO <sub>2</sub> /l
V <sub>talv</sub>	0 to 1000 ml	1 ml for < 100 ml 10 ml for ≥ 100 ml
V <sub>'alv</sub>	0 to 20 l/min	0.01 l/min for < 1 l/min 0.1 l/min for ≥ 1 l/min
V <sub>'CO<sub>2</sub></sub>	0 to 5000 ml/min	1 ml/min for < 1000 ml/min 10 ml/min for ≥ 1000 ml/min
VD <sub>aw</sub>	0 to 1000 ml	1 ml
VD <sub>aw</sub> /VTE	1 to 99 %	1 %
VeCO <sub>2</sub>	0 to 999 ml	1 ml for < 100 ml 10 ml for ≥ 100 ml
ViCO <sub>2</sub>	0 to 999 ml	1 ml for < 100 ml 10 ml for ≥ 100 ml

## Specifications of CO<sub>2</sub> mainstream sensor

**Principle of Operation:** Non-dispersive infrared (NDIR) single beam optics, dual wavelength, no moving parts.

**CO<sub>2</sub> Accuracy:** 0 – 40 mm Hg ± 2 mm Hg  
 41 – 70 mm Hg ± 5% of reading

**Compensations:** Barometric Pressure 400 to 850 mm Hg  
 O<sub>2</sub> and He compensation

**Calibration:** No routine user calibration required. An airway adapter zero is required when changing to a different style of airway adapter.

**Airway Adapters:** < 5 cc deadspace (adult),  
 < 1 cc deadspace (infant)



The blue alarm lamp on top of ventilator lights to indicate heliox is being administered.

## Further reading

Folkow B, Pappenheimer JR. Components of respiratory dead space and their variations with pressure breathing and bronchoactive drugs. *J Appl Physiol* 1955; 8:102-110.

Noe FE. Computer analysis of curves from an infrared CO<sub>2</sub> analyzer and screen-type airflow meter. *J Appl Physiol* 1963; 18:149-157.

Severinghaus JW, Stupfel M. Alveolar dead space as an index of distribution of blood flow in pulmonary capillaries. *J Appl Physiol* 1957; 10:335-348.

Wolff G, Brunner JX, Weibel W. et al. Anatomical and series dead space volume: concept and measurement in clinical practice. *Appl Cardiopul Pathophysiol* 1989; 2:299-307.

Nuckton TJ, Alonso JA, Kallet RH, Daniel BM, Pittet JF, Eisner MD, Matthay MA. Pulmonary dead-space fraction as a risk factor for death in the acute respiratory distress syndrome. *N Engl J Med*. 2002 Apr 25; 346(17):1281-1286.

Astrom E, Niklason L, Drefeldt B, Bajc M, Jonson B. Partitioning of dead space – a method and reference values in the awake human. *Eur Respir J*. 2000 Oct; 16(4):659-664.

